EN: User Manual
 NO: Brukerveiledning
 SV: Bruksanvisning
 ES: Manual del usuario
 PT: Manual do utilizador

Head office:
 GRANBERG AS
 Bjoavegen 1442, NO-5584 Bjoa, Norway
 Tel +47 53 77 53 00, e-mail post@granberg.no

Swedish office / European warehouse:
 GRANBERG SVERIGE AB
 Schubergvägen 20, 311 74 Falkenberg, Sweden
 Tel +46 (0)346 124 25, e-mail post@granberg-ab.se

granberggloves.com

ENGLISH

ART. 210.0093Y – GRANBERG DISPOSABLE ISOLATION GOWN

Description

Disposable isolation gown is made of non-woven polypropylene and polyethylene film. The isolation gowns are 100% latex free. One-size (Universal).

The isolation gowns are made of following material:
 Main material: Polypropylene and polyethylene film.
 Cuffs: Polyester. Sewing thread: Polyester, Polyamide (Nylon).

Dimensions. Length: 132 cm, Chest: 70 cm; Sleeve length: 61 cm

Guide to Apparel and Textile Care Symbols:

EU Type-Examination:

I. This product is classed as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425. This product meets the requirements of the European Standards EN ISO 13688:2013 General Requirements and EN 14126:2003/AC:2004: Protective Clothing - Performance requirements and tests methods for protective clothing against infective agents (Type PB6-B).

PPE CAT. III: Complex design PPE that protects against the risks that may cause very serious consequences such as death or irreversible damage to health.
 II. Notified Body responsible for certification and Module B compliance and ongoing conformity according Module C2 is Autex, Plaza Emilio Sala 1, 03801 ALCOY (ALICANTE), Spain. Notified Body No: 0161.

Performance

Gown performance quoted below is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, mechanical wear, etc.

Levels of protection in compliance with EN 14126:2003/AC:2004:

Properties	Test Method	Results	EN Class*
Abrasion Resistance	EN ISO 12947-2	> 2000 cycles	6 of 6
Flex Cracking Resistance	EN ISO 7854 (Method B)	> 50000 cycles	6 of 6
Trapezoidal Tear Resistance	EN ISO 9073-4	> 20 N	2 of 6
Puncture Resistance	EN 863	> 5 N	1 of 6
Tensile strength	EN ISO 13934-1	> 60N	2 of 6
Resistance to Flame	EN 13274-4 (Method 3)	Pass	1 of 3
Resistance to penetration by synthetic blood under hydrostatic pressure	ISO 16603 (Procedure C)	> 20 kPa	6 of 6
Resistance to penetration by blood borne pathogens using Phi-X174	ISO 16604 (Procedure C)	> 20 kPa	6 of 6
Resistance to Wet Bacterial Penetration**	EN ISO 22610	15< t ≤30	2 of 6
Seam Strength	EN ISO 13935-2	> 50 N	2 of 6
Innocuousness requirements in compliance with the EN ISO 13688:2013			
pH Value	EN ISO 3071	Pass	
AZO Dyes	EN 14362-1	Pass	

* Requirements according to the Standard EN 14325:2018.
 **Information from note * not relevant.

User Guide

Donning Instructions

1. Always remember hand hygiene! Clean your hands before taking the isolation gown on.
2. Start by pulling the individual gown from the bag.
3. Inspect the gown prior to donning for any rips, or holes. If the gown is in good condition, continue to the next step of donning the gown. If there is a hole or tear, discard the gown and restart the donning process.
4. Take hold of the gown from the folded edge and allow the remainder of the gown to unfold. With both hands holding either side of the neck opening swiftly, snap the gown down to release the rest of the gown. Gown is now ready to be donned.
5. Place arms into the sleeves and hands into the thumb-loop cuffs if they apply to your gown.
6. Pull the gown overhead.
7. Pull the side ties to the back of the gown and tie at the rear of the waist.

Doffing Instructions

1. Leaving your medical gloves on, start by crossing your arms grasp on the shoulders and pull the gown away from the body to break the neck closure.
2. Move down to your waist grasp the gown, and pull away from the body to break waist closure.
3. Roll the gown down the body inside out into a ball with your medical gloves to ensure fluids and contaminants are contained.
4. Discard of the gown in an appropriate receptacle. Some users may require special disposal depending on your local facilities rules and applicable regulations.

Packaging Storage & Transportation

Store in a clean room with a relative humidity of not more than 80% and well-ventilated environment with no corrosive gases.
 Preferably stored in original packaging under dry conditions.

The transportation process should be performed in accordance with the contract's provisions. Avoid heavy pressure deformation, direct sunlight, hard objects scratching and rain or snow or liquid leaching.

When stored as recommended, the product will not suffer any chemical properties for up to five years from the date of manufacture. Service life cannot be specified and depends on application and responsibility of user to ascertain suitability of the isolation gown for its intended use.

Use limitations and Disposal:

- Components used in manufacturing may cause allergic reactions in some users. If any allergic reaction occurs, seek medical advice or visit a doctor immediately.
- EU Declaration of Conformity is accessible at www.granberggloves.no/search
- For single use only, Do not re-use.
- The gown does not allow washing.
- Do not use in contact with any types of flame, sparks or heat sources.
- Before usage, inspect isolation gown for any defect or imperfections.
- Handle and disposable the gown in accordance with national regulation.
- The manufacturer cannot be held responsible for improper use of the isolation gown.

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